

Product Datasheet

Recombinant Mouse IgG1 Fc Protein (orb1550467)

Description	Recombinant Mouse IgG1 Fc Protein
Reactivity	Mouse
Endotoxins	< 1 EU/ug by LAL method.
Target	IgG1 Fc
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP2022, Recombinant Mouse IgG1 Fc Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Val98-Lys324) of Mouse IgG1 Fc fused with no tag .
Expiration Date	6 months from date of receipt.