

Product Datasheet

Recombinant Human IL-7R alpha /CD127 Protein, hFc His Tag (orb1550478)

Description	Recombinant Human IL-7R alpha /CD127 Protein, hFc His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	IL7R alpha
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP2011, Recombinant Human IL-7R alpha /CD127 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Glu21-Gly236) of Human IL-7R alpha /CD127 fused with a hFc His Tag tag at the C-terminal.
Expiration Date	6 months from date of receipt.