

Product Datasheet

Recombinant Human FGFR4 Protein, Avi His Tag (orb1550534)

Description	Recombinant Human FGFR4 Protein, Avi His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	FGFR4
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP1955, Recombinant Human FGFR4 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Leu22-Asp369) of Human FGFR4 fused with a His Tag and Avi tag at the C-terminal.
Expiration Date	6 months from date of receipt.