

Product Datasheet

Recombinant Human Adiponectin Protein, His Tag (orb1550598)

Description	Recombinant Human Adiponectin Protein, His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	Adiponectin
Note	For research use only
Purity	> 95% by Tris-Bis PAGE
Source	KMP1891, Recombinant Human Adiponectin Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Glu19-Asn244) of Human Adiponectin fused with a His Tag at the C-terminal.
Expiration Date	6 months from date of receipt.