

Product Datasheet

Recombinant Cyno Cyno ACE2 Protein, Avi His Tag (orb1550600)

Description	Recombinant Cyno Cyno ACE2 Protein, Avi His Tag
Endotoxins	< 1 EU/ug by LAL method.
Target	Cyno Cyno ACE2
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP1889, Recombinant Cyno Cyno ACE2 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Gln18-Val739) of Cyno Cyno ACE2 fused with a His Tag and Avi at the C-terminal.
Expiration Date	6 months from date of receipt.