

Product Datasheet

Recombinant Streptococcus pyogenes M1 CRISPR-Cas9 Protein, His Tag (orb1550691)

Description	Recombinant Streptococcus pyogenes M1 CRISPR-Cas9 Protein, His Tag
Species/Host	Virus
Reactivity	Bacteria
Endotoxins	< 1.0 EU per ug proteinby the LAL method.
Target	M1 CRISPRCas9
Note	For research use only
Purity	> 90 %by SDS-PAGE
Source	KMP1798, Recombinant Streptococcus pyogenes M1 CRISPR-Cas9 Protein is produced by Baculovirus-Insect Cells expression system. The target protein is expressed with sequence (Met1-Asp1368) of Streptococcus pyogenes M1 CRISPR-Cas9 (Accession #WP_010922251.1) fused with a nuclear localization sequence (NLS) at the N-terminus, two NLSs and a polyHis Tagtidine tag at the C-terminus.
Expiration Date	6 months from date of receipt.