

## Product Datasheet

# Recombinant Human ACE2 Protein with hFc tag, Fc His Tag (orb1550698)

|                        |                                                                                                                                                                                                                              |
|------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | Recombinant Human ACE2 Protein with hFc tag, Fc His Tag                                                                                                                                                                      |
| <b>Species/Host</b>    | Human                                                                                                                                                                                                                        |
| <b>Reactivity</b>      | Human                                                                                                                                                                                                                        |
| <b>Endotoxins</b>      | < 0.1 EU/ug by LAL method.                                                                                                                                                                                                   |
| <b>Target</b>          | ACE2                                                                                                                                                                                                                         |
| <b>Note</b>            | For research use only                                                                                                                                                                                                        |
| <b>Purity</b>          | > 95% by SDS-PAGE.                                                                                                                                                                                                           |
| <b>Source</b>          | KMP1791, Recombinant Human ACE2 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Gln18-Ser740) of human ACE2 (Accession #Q9BYF1) fused with a Fc tag at the C-terminus. |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                                                                                                                                                               |