

Product Datasheet

Recombinant human FGL1 Protein, His Tag (orb1550699)

Description	Recombinant human FGL1 Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	Please contact us for more information.
Target	FGL1
Note	For research use only
Purity	> 90% by SDS-PAGE.
Source	KMP1790, Recombinant human FGL1 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Met1-Ile312) of human FGL1 (Accession #NP_004458.3) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.