

Product Datasheet

Recombinant Human TPBG Protein, His Tag (orb1550717)

Description	Recombinant Human TPBG Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	TPBG
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1772, Recombinant Human TPBG Protein is produced by Mammalian expression system. The target protein is expressed with sequence (Ser32-Ser355) of human TPBG (Accession #Q13641) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.