

Product Datasheet

Recombinant Human Semaphorin-5A Protein, His Tag (orb1550719)

Description	Recombinant Human Semaphorin-5A Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	Semaphorin5A
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1770, Recombinant Human Semaphorin-5A Protein is produced by Mammalian expression system. The target protein is expressed with sequence (Glu23-Thr765) of human Semaphorin-5A (Accession #Q13591) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.