

Product Datasheet

Recombinant Human AZGP1 Protein, His Tag (orb1550722)

Description	Recombinant Human AZGP1 Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	AZGP1
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1767, Recombinant Human AZGP1 Protein is produced by Mammalian expression system. The target protein is expressed with sequence (Gln21-Ser298) of human AZGP1 (Accession #P25311) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.