

Product Datasheet

Recombinant Human Adiponectin/Acrp30/ADIPOQ Protein, His Tag (orb1550731)

Description	Recombinant Human Adiponectin/Acrp30/ADIPOQ Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	Adiponectin
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1758, Recombinant Human Adiponectin/Acrp30/ADIPOQ Protein is produced by Human Cells expression system. The target protein is expressed with sequence (Glu19-Asn244) of human Adiponectin/Acrp30/ADIPOQ (Accession #Q15848) fused with a 6xHis Tag at the C-terminus.
Expiration Date	6 months from date of receipt.