

Product Datasheet

Recombinant Mouse VEGF 164 Protein, His Tag (orb1550780)

Description	Recombinant Mouse VEGF 164 Protein, His Tag
Species/Host	Human
Reactivity	Mouse
Endotoxins	< 0.1 EU/ug by LAL method.
Target	VEGF 164
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1709, Recombinant Mouse VEGF 164 Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Ala27-Arg190) of mouse VEGF 164 (Accession #NP_001273986.1.) fused with a 6×His Tag at the N-terminus.
Expiration Date	6 months from date of receipt.