

Product Datasheet

Recombinant Human CCL18/PARC Protein, His Tag (orb1550816)

Description	Recombinant Human CCL18/PARC Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	CCL18
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1673, Recombinant Human CCL18/PARC Protein is produced by E. coli expression system. The target protein is expressed with sequence (Ala21-Ala89) of human CCL18/PARC (Accession #P55774) fused with an initial Met at the N-terminus and a 6×His Tag at the N-terminus.
Expiration Date	6 months from date of receipt.