

Product Datasheet

Recombinant Rat PD-L1/B7-H1/CD274 Protein, His Tag (orb1550817)

Description	Recombinant Rat PD-L1/B7-H1/CD274 Protein, His Tag
Species/Host	Human
Reactivity	Rat
Endotoxins	< 1 EU/ug by LAL method.
Target	PDL1
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1672, Recombinant Rat PD-L1/B7-H1/CD274 Protein is produced by Human Cells expression system. The target protein is expressed with sequence (Ala18-Thr238) of rat PD-L1/B7-H1/CD274 (Accession #D4AE25) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.