

Product Datasheet

Recombinant Mouse IL-15RA/IL-15 R alpha Protein, Fc His Tag (orb1550830)

Description	Recombinant Mouse IL-15RA/IL-15 R alpha Protein, Fc His Tag
Species/Host	Human
Reactivity	Mouse
Endotoxins	< 1 EU/ug by LAL method.
Target	IL15RA
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1659, Recombinant Mouse IL-15RA/IL-15 R alpha Protein is produced by Human Cells expression system. The target protein is expressed with sequence (Gly33-Lys205) of mouse IL-15RA/IL-15 R alpha (Accession #Q60819) fused with an Fc tag at the C-terminus.
Expiration Date	6 months from date of receipt.