

Product Datasheet

Recombinant Human FABP3/H-FABP Protein, His Tag (orb1550862)

Description	Recombinant Human FABP3/H-FABP Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	FABP3
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1627, Recombinant Human FABP3/H-FABP Protein is produced by E. coli expression system. The target protein is expressed with sequence (Val2-Ala133) of human FABP3/H-FABP (Accession #P05413) fused with an initial Met at the N-terminus and a 6×His Tag at the N-terminus.
Expiration Date	6 months from date of receipt.