

Product Datasheet

Recombinant Human TrkB Protein, His Tag (orb1550896)

Description	Recombinant Human TrkB Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	TrkB
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1593, Recombinant Human TrkB Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Cys32-His Tag430) of human TrkB (Accession #NP_001018074.1) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.