

Product Datasheet

Recombinant Human ICAM-2 Protein, Fc & His Tag (orb1550913)

Description	Recombinant Human ICAM-2 Protein, Fc & His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	Please contact us for more information.
Target	ICAM2
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1576, Recombinant Human ICAM-2 Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Lys25-Gln216) of human ICAM-2 (Accession #NP_000864.2) fused with an Fc, 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.