

Product Datasheet

Recombinant Human CHIT1 Protein, His Tag (orb1550950)

Description	Recombinant Human CHIT1 Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.01 EU/ug by LAL method.
Target	CHIT1
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1538, Recombinant Human CHIT1 Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Ala22-Asn466) of human CHIT1 (Accession #NP_003456) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.