

Product Datasheet

Recombinant Human L-Selectin/CD62L Protein, Fc & His Tag (orb1550953)

Description	Recombinant Human L-Selectin/CD62L Protein, Fc & His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	LSelectin
Note	For research use only
Purity	> 92% by SDS-PAGE.
Source	KMP1535, Recombinant Human L-Selectin/CD62L Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Trp52-Asn345) of human L-Selectin/CD62L (Accession #NP_000646.2) fused with an Fc, 6xHis Tag at the C-terminus.
Expiration Date	6 months from date of receipt.