

Product Datasheet

Recombinant Human CCL17/TARC Protein, His Tag (orb1550996)

Description	Recombinant Human CCL17/TARC Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	CCL17
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1491, Recombinant Human CCL17/TARC Protein is produced by Human Cell expression system. The target protein is expressed with sequence (Ala24Ser94) of human CCL17/TARC (Accession #Q92583) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.