

Product Datasheet

Recombinant Human LAIR2 Protein, His Tag (orb1551064)

Description	Recombinant Human LAIR2 Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	LAIR2
Note	For research use only
Purity	> 90% by SDS-PAGE.
Source	KMP1423, Recombinant Human LAIR2 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Gln22-Pro152) of human LAIR2 (Accession #NP_002279) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.