

Product Datasheet

Recombinant Human CEACAM6 Protein, His Tag (orb1551071)

Description	Recombinant Human CEACAM6 Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	CEACAM6
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1416, Recombinant Human CEACAM6 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Lys35-Gly320) of human CEACAM6 (Accession #NP_002474.3) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.