

Product Datasheet

Recombinant Human ENPP-2 Protein, His Tag (orb1551223)

Description	Recombinant Human ENPP-2 Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	ENPP2
Note	For research use only
Purity	> 90% by SDS-PAGE.
Source	KMP1264, Recombinant Human ENPP-2 Protein is produced by Mammalian expression system. The target protein is expressed with sequence (Ala36-Ile863) of human ENPP-2 (Accession #AAH34961.1) fused with a 6xHis Tag at the C- terminus.
Expiration Date	6 months from date of receipt.