

Product Datasheet

Recombinant Human ACVR2B/Activin RIIIB Protein, Fc & His Tag (orb1551312)

Description	Recombinant Human ACVR2B/Activin RIIIB Protein, Fc & His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	ACVR2B
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1175, Recombinant Human ACVR2B/Activin RIIIB Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Ser19-Thr134) of human ACVR2B/Activin RIIIB (Accession #NP_001097.2) fused with an Fc, 6xHis Tag at the C-terminus.
Expiration Date	6 months from date of receipt.