

Product Datasheet

Recombinant Human Carbonic Anhydrase XII/CA12 Protein, His Tag (orb1551334)

Description	Recombinant Human Carbonic Anhydrase XII/CA12 Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	Carbonic Anhydrase XII
Note	For research use only
Purity	> 90% by SDS-PAGE.
Source	KMP1153, Recombinant Human Carbonic Anhydrase XII/CA12 Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Ala25-Gln291) of human Carbonic Anhydrase XII/CA12 (Accession #NP_001209.1) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.