

Product Datasheet

Recombinant Human CD89/FCAR Protein, Fc & His Tag (orb1551365)

Description	Recombinant Human CD89/FCAR Protein, Fc & His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	CD89
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1122, Recombinant Human CD89/FCAR Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Gln22-Asn227) of human CD89/FCAR (Accession #NP_001991.1) fused with an Fc, 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.