

Product Datasheet

Recombinant Human IL-1 beta Protein, His Tag (orb1551467)

Description	Recombinant Human IL-1 beta Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	IL1 beta
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1020, Recombinant Human IL-1 beta Protein is produced by E. coli expression system. The target protein is expressed with sequence (Ala117-Ser269) of human IL-1 beta (Accession #NP_000567.1) fused with an initial Met at the N-terminus and a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.