

Product Datasheet

Recombinant Human TFIIB/GTF2B Protein, Gst Tag (orb1551563)

Description	Recombinant Human TFIIB/GTF2B Protein, Gst Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	GTF2B
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human TFIIB(Met1-Leu316) was fused with the N-terminal Gst Tag and expressed in E. coli.
Uniprot ID	Q00403
Expiration Date	6 months from date of receipt.