

Product Datasheet

Recombinant Human TFF2 Protein, His Tag (orb1551564)

Description	Recombinant Human TFF2 Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	TFF2
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human TFF2(Gln24-Tyr129) was fused with the C-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	Q03403
Expiration Date	6 months from date of receipt.