

Product Datasheet

Recombinant Mouse Testican-3/Spock3 Protein, His Tag (orb1551566)

Description	Recombinant Mouse Testican-3/Spock3 Protein, His Tag
Species/Host	E. coli
Reactivity	Mouse
Endotoxins	Please contact our customer service for more information.
Target	Spock3
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the mouse Spock3(Ala23-Ile436) was fused with the C-terminal His Tag and expressed in E. coli.
Uniprot ID	Q8BKV0
Expiration Date	6 months from date of receipt.