

Product Datasheet

Recombinant Human Syntaxin-7/STX7 Protein, His Tag (orb1551575)

Description	Recombinant Human Syntaxin-7/STX7 Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	STX7
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human STX7(Ser2-Leu238) was fused with the N-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	O15400
Expiration Date	6 months from date of receipt.