

Product Datasheet

Recombinant Human SUMO3 AMC Protein, No Tag (orb1551582)

Description	Recombinant Human SUMO3 AMC Protein, No Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	SUMO3
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human SUMO3(Met1-Phe103) was fused without Tag and expressed in Mammalian cells.
Uniprot ID	P55854
Expiration Date	6 months from date of receipt.