

Product Datasheet

Recombinant Human SUMO1 Protein, His Tag (orb1551584)

Description	Recombinant Human SUMO1 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	SUMO1
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human SUMO1(Met1-Val101) was fused with the N-terminal His Tag and expressed in E. coli.
Uniprot ID	P63165
Expiration Date	6 months from date of receipt.