

Product Datasheet

Recombinant Human SULT1C4 Protein, His Tag (orb1551589)

Description	Recombinant Human SULT1C4 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	SULT1C4
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human SULT1C4(Met1-Lys102) was fused with the N-terminal His Tag and expressed in E. coli.
Uniprot ID	075897
Expiration Date	6 months from date of receipt.