

## Product Datasheet

### Recombinant Human SERPINA1 Protein, His Tag (orb1551686)

|                        |                                                                                                                                |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | Recombinant Human SERPINA1 Protein, His Tag                                                                                    |
| <b>Species/Host</b>    | Mammal                                                                                                                         |
| <b>Reactivity</b>      | Human                                                                                                                          |
| <b>Endotoxins</b>      | < 1.0 EU per µg of the protein as determined by the LAL method.                                                                |
| <b>Target</b>          | SERPINA1                                                                                                                       |
| <b>Note</b>            | For research use only                                                                                                          |
| <b>Purity</b>          | ≥90% as determined by SDS-PAGE                                                                                                 |
| <b>Source</b>          | This protein contains the human SERPINA1(Glu25-Lys418) was fused with the C-terminal His Tag and expressed in Mammalian cells. |
| <b>Uniprot ID</b>      | <b>P01009</b>                                                                                                                  |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                                                                 |