

Product Datasheet

Recombinant Human SENP2 Protein, His & SUMO Tag (orb1551690)

Description	Recombinant Human SENP2 Protein, His & SUMO Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	SENP2
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human SENP2(Asp363-Leu589) was fused with the N-terminal His and SUMO Tag and expressed in E. coli.
Uniprot ID	Q9HC62
Expiration Date	6 months from date of receipt.