

Product Datasheet

Recombinant Human Semaphorin 5A/SEMA5A Protein, His Tag (orb1551692)

Description	Recombinant Human Semaphorin 5A/SEMA5A Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	SEMA5A
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human SEMA5A(Glu23-Thr765) was fused with the C-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	Q13591
Expiration Date	6 months from date of receipt.