

Product Datasheet

EDRF1 antibody (orb1552)

Description	EDRF1 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthetic peptide derived from Human EDRF1/C10orf137 (101-200/1238 aa)
Target	EDRF1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	139 kDa
Uniprot ID	Q3B7T1
Expiration Date	12 months from date of receipt.