

Product Datasheet

Recombinant Human MAPKSP1/LAMTOR3 Protein, His Tag (orb1552047)

Description	Recombinant Human MAPKSP1/LAMTOR3 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	LAMTOR3
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human LAMTOR3(Met1-Val124) was fused with the N-terminal His Tag and expressed in E. coli.
Uniprot ID	Q9UHA4
Expiration Date	6 months from date of receipt.