

Product Datasheet

Recombinant Mouse Light/Tnfsf14/CD258 Protein, His Tag (orb1552090)

Description	Recombinant Mouse Light/Tnfsf14/CD258 Protein, His Tag
Species/Host	Mammal
Reactivity	Mouse
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	Tnfsf14
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the mouse Tnfsf14(Asp72-Val239) was fused with the N-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	Q9QYH9
Expiration Date	6 months from date of receipt.