

## Product Datasheet

### Recombinant Human LGALSL Protein, His Tag (orb1552094)

|                        |                                                                                                                              |
|------------------------|------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | Recombinant Human LGALSL Protein,<br>His Tag                                                                                 |
| <b>Species/Host</b>    | E. coli                                                                                                                      |
| <b>Reactivity</b>      | Human                                                                                                                        |
| <b>Endotoxins</b>      | Please contact our customer service<br>for more information.                                                                 |
| <b>Target</b>          | LGALSL                                                                                                                       |
| <b>Note</b>            | For research use only                                                                                                        |
| <b>Purity</b>          | ≥90% as determined by SDS-PAGE                                                                                               |
| <b>Source</b>          | This protein contains the human<br>LGALSL(Ala2-Ser172) was fused with<br>the C-terminal His Tag and expressed<br>in E. coli. |
| <b>Uniprot ID</b>      | <b>Q3ZCW2</b>                                                                                                                |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                                                               |