

Product Datasheet

Recombinant Mouse IL-18/IL-1F4 Protein, His Tag (orb1552184)

Description	Recombinant Mouse IL-18/IL-1F4 Protein, His Tag
Species/Host	Mammal
Reactivity	Mouse
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	Il18
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the mouse Il18(Asn36-Ser192) was fused with the C-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	P70380
Expiration Date	6 months from date of receipt.