

Product Datasheet

Recombinant Human EpCAM/TROP1 Protein, Fc Tag (orb1552493)

Description	Recombinant Human EpCAM/TROP1 Protein, Fc Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	EPCAM
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human EPCAM(Gln24-Lys265) was fused with the C-terminal Fc Tag and expressed in Mammalian cells.
Uniprot ID	P16422
Expiration Date	6 months from date of receipt.