

Product Datasheet

Recombinant Human DTYMK Protein, His Tag (orb1552523)

Description	Recombinant Human DTYMK Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	DTYMK
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human DTYMK(Met1-Lys212) was fused with the C-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	P23919
Expiration Date	6 months from date of receipt.