

## Product Datasheet

### Recombinant Human CXCL11/I-TAC Protein, No Tag (orb1552592)

|                        |                                                                                                     |
|------------------------|-----------------------------------------------------------------------------------------------------|
| <b>Description</b>     | Recombinant Human CXCL11/I-TAC Protein, No Tag                                                      |
| <b>Species/Host</b>    | E. coli                                                                                             |
| <b>Reactivity</b>      | Human                                                                                               |
| <b>Endotoxins</b>      | Please contact our customer service for more information.                                           |
| <b>Target</b>          | CXCL11                                                                                              |
| <b>Note</b>            | For research use only                                                                               |
| <b>Purity</b>          | ≥90% as determined by SDS-PAGE                                                                      |
| <b>Source</b>          | This protein contains the human CXCL11(Phe22-Phe94) was fused without Tag and expressed in E. coli. |
| <b>Uniprot ID</b>      | <b>O14625</b>                                                                                       |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                                      |