

Product Datasheet

Recombinant Human CXCL10/IP-10 Protein, His & Fc Tag (orb1552594)

Description	Recombinant Human CXCL10/IP-10 Protein, His & Fc Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	CXCL10
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human CXCL10(Val22-Pro98) was fused with the C-terminal His Tag and Fc Tag and expressed in Mammalian cells.
Uniprot ID	P02778
Expiration Date	6 months from date of receipt.