

## Product Datasheet

### Recombinant Mouse Cx3cl1 Protein, His Tag (orb1552597)

<b>Description</b>	Recombinant Mouse Cx3cl1 Protein, His Tag
<b>Species/Host</b>	Mammal
<b>Reactivity</b>	Mouse
<b>Endotoxins</b>	< 1.0 EU per µg of the protein as determined by the LAL method.
<b>Target</b>	Cx3cl1
<b>Note</b>	For research use only
<b>Purity</b>	≥90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the mouse Cx3cl1(Gln25-Arg337) was fused with the C-terminal His Tag and expressed in Mammalian cells.
<b>Uniprot ID</b>	<b>O35188</b>
<b>Expiration Date</b>	6 months from date of receipt.