

Product Datasheet

Recombinant Human CKMM/CKM Protein, His Tag (orb1552678)

Description	Recombinant Human CKMM/CKM Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 1.0 EU per µg of the protein as determined by the LAL method.
Target	CKM
Note	For research use only
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human CKM(Met1-Lys381) was fused with the C-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	P06732
Expiration Date	6 months from date of receipt.